



Intense THz generation with 1300-1500 nm light

DAST 100-fs pump, 0.84 mJ 1450-nm pump





High-quality DAST has been the work-horse organic crystal for broadband, high-field THz genereation for several years. With a range of crystal thicknesses and extremely large aperture sizes available, we can meet your high-field THz generation needs. DAST is optimally phase-matched for pumping from 1.3-1.6 microns, but also provides good phase-matching from 1-2 microns.

DAST vs OH1 vs BNA-S

100-fs pump

BNA-S

4.2 mJ 800-nm pump vs. DAST, OH1

o.84 mJ 1450-nm pump (same THz focusing setup)

